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Instruction for use

1. Product Description

optimill memosplint and **memosplint Version S** products are resin discs for the CAD/CAM process for the production of flexible custom-made occlusal splints by dental specialists. Due to its greater hardness, **optimill memosplint Version S** can also be used for interim restorations with clamps. The milled works have a thermoelastic flexibility.
Patient target group: adults and adolescents

2. Intended purpose

optimill memosplint:

Dental resins for the production of flexible dental therapeutic occlusal splints of all kinds.

optimill memosplint Version S:

Dental resins for the production of flexible dental therapeutic occlusal splints of all kinds and interim restorations with clamps.

3. Contraindications

The material should not be used for any purpose other than the generative CAD/CAM fabrication of dental bite splints. The material is not suitable for fabricating total denture bases. Undesirable biological reactions (e.g. allergies to material components) may occur in very rare cases. Do not use bite splints made with this material if incompatibilities or allergies to the material components are known.

4. Risk minimization and safety instructions


Processing the disc produces dust that may irritate the eyes, skin and respiratory tract. So please always ensure that the extractor fan at your workplace is working flawlessly. See safety data sheet on www.dentona.de.

5. Processing steps

5.1 Designing

When splints are being designed, the wall thickness (also of the clamps for **optimill memosplint Version S**) must not be less than the minimum of 1.0 mm. Please check this requirement in your design software. Any corners and edges must be rounded off.

During nesting, please ensure that enough attachments are used. It is a thermoelastic material and would otherwise become deformed in the machine. This would result in tool breakage.

Nesting **optimill memosplint**

Nesting **optimill memosplint Version S**

5.2 Milling

A conventional dry milling strategy for PMMA can be used. Single-edged PMMA tools should be used in the machine. The attachments are cleaned as usual after milling.

5.3 Finishing

The works can be processed with suitable cross-toothed carbide cutters or suitable cutting wheels for resin. Final polishing of the splints is essential to prevent any accumulation of plaque on them. Polishing is done with standard polishing brushes and polishing pastes. Attention: Excessive contact pressure during polishing may overheat and melt the material.

6. Cleaning

The objects should not be cleaned with excessive evaporation due to the high heat generation. Clean the objects for approx. 5 minutes in water in an ultrasonic bath at a temperature of 35°C until no milling or polishing residue can be seen or felt. Cleaning solutions, alcoholic disinfectants and denture cleaners are not recommended.

7. Notes for the laboratory / for the dentist

Restorations are custom-made products according to EU regulations. Store in clear water overnight before delivery. For each process, indicate the LOT number required to identify the material. Communicate the below handling information to the patient.

7.1 Notes for the patient

For better comfort, the splint can be placed in 30 - 40 °C warm water just before insertion. After wearing, clean under running water with a soft toothbrush without toothpaste. Do not use prosthesis cleaners, cleaning solutions or alcohol-based disinfectants for cleaning. Then keep it dry until next wear. Protect the splint from excessive heat and cold and avoid direct sunlight.

8. Technical Data

Parametres	optimill memosplint	optimill memosplint Version S
Shore D (21° – 25°C)	> 60	> 65
Shore D (36° – 38° C)	< 60	< 65
Density	approx. 1.1-1.2 g/cm ³	
Color	clear transparent	

9. Chemical composition in %

Poly(m)ethylacrylate - cross-linked copolymers of methacrylic acid > 90

10. Storage

Avoid direct sunlight. No further special measures required.

Do not use after the expiration date.

11. Disposal

Restorations and remains of this material are insoluble in water, inactive, pose no threat to ground water and can thus be disposed of in normal domestic waste.

12. Information

Our application-specific recommendations, no matter whether they are given orally, in writing or by way of practical tutorials, are based on our own experiences and tests and may therefore only be considered as recommended values. Our products are subject to continuous enhancements. We therefore reserve the right to make changes in design and composition.

Report all serious incidents (death, serious deterioration of health, serious risk to public health) that have occurred in connection with the medical device to the manufacturer and the competent authority of the member state.

13. Symbols

The symbols used in the labeling are explained on the last page.

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dentona AG | Otto-Hahn-Str. 27 | DE-44227 Dortmund